



New Dimethachlor metabolites


The European Food Safety Authority (EFSA) listed in the „Scientific Report 2008“, the relevant metabolites of the pesticide Dimethachlor. Beside the two metabolites Dimethachlor OA and ESA we can now offer you the following new metabolites: CGA39981, CGA42443, CGA102935, CGA369873 and SYN528702. All metabolites have a purity >95% and will be delivered with a certificate of analysis.




Product name	Dimethachlor metabolite CGA 39981
Synonyme:	N-(2,6-dimethylphenyl)-2-hydroxy-N-(2-methoxyethyl)acetamide
Article No.	RCG-311
CAS-No	-
Formula	C ₁₃ H ₁₉ NO ₃
Molecular weight	237.30g/mol
Purity	≥95%
Pack size:	100mg à 310CHF
	Bulk quantities on request




Product name	Dimethachlor metabolite CGA 42443
Synonyme:	N-(2,6-dimethylphenyl)-N-(2-methoxyethyl)acetamide
Article No.	RCG-312
CAS-No	-
Formula	C ₁₃ H ₁₉ NO ₂
Molecular weight	221.30g/mol
Purity	≥95%
Pack size:	100mg à 250CHF
	Bulk quantities on request




Product name	Dimethachlor metabolite CGA 102935
Synonyme:	N-carboxymethyl-N-(2,6-dimethylphenyl)oxalamic acid
Article No.	RCG-313
CAS-No	-
Formula	C ₁₂ H ₁₃ NO ₅
Molecular weight	251.24g/mol
Purity	≥95%
Pack size:	100mg à 310CHF
	Bulk quantities on request






Product name	Dimethachlor metabolite CGA-369873
Synonyme:	(2,6-dimethylphenylcarbamoyl)-methanesulfonic acid sodium salt
Article No.	RCG-314
CAS-No	-
Lot No	2009441L4.4
Formula	C ₁₀ H ₁₂ NNaO ₄ S
Molecular weight	256.27 g/mol
Purity	≥95%
Pack size:	100mg à 210CHF
	Bulk quantities on request



Product name	Dimethachlor OA
CAS-No	-
Article No.	RCG-208
Formula	C ₁₃ H ₁₇ NO ₄
Molecular weight	251.28 g/mol
Purity	≥95%
Pack size:	10mg à 180CHF
	Bulk quantities on request



Product name	Dimethachlor ESA sodium salt
CAS-No	-
Article No.	RCG-209
Formula	C ₁₃ H ₁₈ NNaO ₅ S
Molecular weight	323.35g/mol
Purity	≥95%
Pack size:	10mg à 180CHF
	Bulk quantities on request

If you are interested in other metabolites, please do not hesitate to contact us. With our broad knowledge in the synthesis of metabolites and other reference materials we are the reliable partner for you.

We help you make your molecules!