





New Reference standards:

- 4-Methylmethcathinone hydrochloride
- 4-Methoxymethcathinone hydrochloride
- Butylone hydrochloride


Substituted cathinones, like 4-Methylmethcathinone and 4-Methoxymethcathinone, and also the β -keto analogue of MBDB (Butylone hydrochloride) are available in various online shops as research chemicals. ReseaChem has synthesized a first set of those „research chemicals“ as reference standards. The standards are available in deuterated and non-deuterated form.



| | |
|------------------|--|
| Product name | 4-Methylmethcathinone hydrochloride |
| Synonyme: | Mephedrone, 4-MMC, 4-Methylephedrone |
| Article No. | RCG-301 |
| CAS-No | 1189726-22-4 |
| Formula | $C_{11}H_{15}NO \cdot HCl$ |
| Molecular weight | 213.71g/mol |
| Purity | >98.5% |
| Pack size: | 10MG CHF 115.-- RCG-301-10MG Bulk quantities on request |




| | |
|------------------|--|
| Product name | 4-Methylmethcathinone-d ₃ hydrochloride |
| Synonyme: | |
| Article No. | RCG-322 |
| CAS-No | |
| Formula | $C_{11}H_{12}NOD_3 \cdot HCl$ |
| Molecular weight | 216.73g/mol |
| Purity | >98.5% |
| Pack size: | 5MG CHF 120.-- RCG-322-5MG Bulk quantities on request |




| | |
|------------------|--|
| Product name | 4-Methoxymethcathinone hydrochloride |
| Synonyme: | para-Methoxymethcathinone, bk-PMMA, PMMC, Methedrone |
| Article No. | RCG-321 |
| CAS-No | 530-54-1 |
| Formula | $C_{11}H_{15}NO \cdot HCl$ |
| Molecular weight | 229.70g/mol |
| Purity | >98.5% |
| Pack size: | 10MG CHF 115.-- RCG-321-10MG Bulk quantities on request |






Product name 4-Methoxymethcathinone-d₃ hydrochloride
 Synonyme:
 Article No. RCG-323
 CAS-No
 Formula C₁₁H₁₂NO₂D₃ . HCl
 Molecular weight 232.73g/mol
 Purity >98.5%
 Pack size: 5MG CHF 120.-- RCG-323-5MG
 Bulk quantities on request



Product name Butylone hydrochloride
 Synonyme: bk-MBDB
 Article No. RCG-302
 CAS-No 17762-90-2
 Formula C₁₂H₁₅NO₃ . HCl
 Molecular weight 257.72g/mol
 Purity >98.5%
 Pack size: 10MG CHF 115.-- RCG-302-10MG
 Bulk quantities on request



Product name Butylone-d₃ hydrochloride
 Synonyme:
 Article No. RCG-324
 CAS-No
 Formula C₁₂H₁₂NO₃D₃ . HCl
 Molecular weight 260.74g/mol
 Purity >98.5%
 Pack size: 5MG CHF 120.-- RCG-324-5MG
 Bulk quantities on request

Coming soon:

Methylone, the the β-keto analogue of MDMA, deuterated and non-deuterated.

Missing a compound?

If you are interested in other compounds, please do not hesitate to contact us.