


1-(3-Chlorophenyl)-piperazine-D8 hydrochloride (mCPP-D8)

meta-Chlorophenylpiperazine (mCPP) is a psychoactive drug of the phenylpiperazine class. It's often sold as ecstasy (XTC).

We are proud to present you deuterated **1-(3-Chlorophenyl)-piperazine hydrochloride** as reference standard for GC-MS and LC-MS.



Product name	1-(3-Chlorophenyl)-piperazine-D8 hydrochloride		
Article No.	RCG-372		
Formula	C ₁₀ H ₅ ClN ₂ D ₈ . HCl		
Molecular weight	241.19g/mol		
Purity	≥98.5% (HPLC)		
Pack size:	RCG-372-10MG	10MG	CHF 180.--
	RCG-372-100MG	100MG	CHF 1'550.--
	RCG-372-1G	1G	CHF 4'500.--

other quantities on request

If you are interested in other reference standard material, please do not hesitate to contact us. With our broad knowledge in the synthesis of metabolites and other reference materials we are the reliable partner for you.

We help to make your molecules!

