

NMR Service Price List

as of March 2026

ReseaChem offers a unique 24-hour-NMR-Service for ^1H -NMR analysis. Additionally, we offer a variety of other NMR-Services on a brand-new Bruker 400MHz NMR-spectrometer (AvanceCore).

Due to rising costs in materials, maintenance, and improved spectral resolution, ReseaChem is forced to adjust its pricing for the NMR-Service per March 1st 2026 as follows:

1D-NMR

For measurements in common solvents (CDCl_3 , DMSO-D_6 , D_2O , Methanol-D_3 , Acetonitrile-D_3 , Acetone-D_6) we charge the following rates:

- ^1H -NMR (within 24 hours): CHF 130.-
- ^{13}C -NMR: CHF 150.- *
- ^{13}C -DEPT-NMR: CHF 130.-
- ^{13}C -APT-NMR: CHF 130.-
- ^{31}P -NMR: CHF 130.-
- ^{19}F -NMR: CHF 130.-
- Quantitative ^1H -NMR**: CHF 220.-
- Quantitative ^{31}P -NMR**: CHF 220.-
- Quantitative ^{19}F -NMR**: CHF 220.-

*The measurement is performed with 1000 scans (Ideally 1 mg per carbon atom in the molecule is delivered to obtain good spectral quality). Additional acquisition time will be charged (CHF 130/1000 scans).

**If a determination for the internal standard in quantitative NMR is necessary, we additionally charge CHF 130.- for the ^1H -NMR. Otherwise, the customer is responsible for the corresponding choice of an appropriate internal standard.

2D-NMR

For measurements in common solvents (CDCl_3 , DMSO-D_6 , D_2O , Methanol-D_4 , Acetonitrile-D_3 , Acetone-D_6) we charge the following rates:

- $^1\text{H}/^1\text{H}$ -COSY ($^3\text{J}_{\text{H,H}}$): CHF 230.-
- $^{13}\text{C}/^1\text{H}$ -HSQC ($^1\text{J}_{\text{C,H}}$): CHF 230.-
- Multiplicity Edited HSQC: CHF 230.-
- $^{13}\text{C}/^1\text{H}$ -HMBC ($^2,^3\text{J}_{\text{C,H}}$): CHF 230.-
- $^1\text{H}/^1\text{H}$ -NOESY: CHF 230.-
- $^1\text{H}/^1\text{H}$ -ROESY: CHF 230.-

***2D-NMR does not include the separate ^1H and ^{13}C -spectra – if necessary they will be charged accordingly.

Other NMR-Experiments on demand.

Uncommon Solvents (Surcharge per dissolved sample)

- Benzene- D_6 : CHF 30.-
- Methanol- D_4 : CHF 30.-
- Methylene chloride- D_2 : CHF 40.-
- N,N-DMF- D_7 : CHF 50.-
- THF- D_8 : CHF 50.-
- Toluene- D_8 : CHF 30.-

Other solvents on demand.

If you have any questions or special requests, please don't hesitate to contact us.